

Accumax Tips Comappatibility Chart vs Accumax FAB Pipette

Non-Filter tips				Compatible Volumes (μ l)		
Packing style:	Bulk	Rack	Reloading stack			
Volume Range:	10 μ l			Variable	Fixed	Multi
	Catalog #					
Standard	AT-10-B	AT-10-R	AT-10-K	0.1 - 2.5 μ l 0.5-10 μ l	2.5 μ l	0.5 - 10 μ l 2-20 μ l
Sterile	AT-10-S-B	AT-10-S-R	AT-10-S-K			
Sterile, Low retention	AT-10-S-L-B	AT-10-S-L-R	AT-10-S-L-K			
Packing style:	Bulk	Rack	Reloading stack			
Volume Range:	10XL			Variable	Fixed	Multi
	Catalog #					
Standard	AT-10XL-B	AT-10XL-R	AT-10XL-K	0.1 - 2.5 μ l 0.5-10 μ l	2.5 μ l	0.5 - 10 μ l 2-20 μ l
Sterile	AT-10XL-S-B	AT-10XL-S-R	AT-10XL-S-K			
Sterile, Low retention	AT-10XL-S-L-B	AT-10XL-S-L-R	AT-10XL-S-L-K			
Packing style:	Bulk	Rack	Reloading stack			
Volume Range:	10XN			Variable	Fixed	Multi
	Catalog #					
Standard	AT-10XN-B	AT-10XN-R	AT-10XN-K	0.1 - 2.5 μ l 0.5-10 μ l	2.5 μ l	0.5 - 10 μ l 2-20 μ l
Sterile	AT-10XN-S-B	AT-10XN-S-R	AT-10XN-S-K			
Sterile, Low retention	AT-10XN-S-L-B	AT-10XN-S-L-R	AT-10XN-S-L-K			
Packing style:	Bulk	Rack	Reloading stack			
Volume Range:	200 μ l			Variable	Fixed	Multi
	Catalog #					
Standard	AT-200-B	AT-200-R	AT-200-K	2-20 μ l 5-50 μ l 10-100 μ l 20-200 μ l	5 μ l, 10 μ l, 20 μ l, 25 μ l, 50 μ l, 100 μ l, 200 μ l	5- 50 μ l 10-100 μ l 20-200 μ l
Sterile	AT-200-S-B	AT-200-S-R	AT-200-S-K			
Sterile, Low retention	AT-200-S-L-B	AT-200-S-L-R	AT-200-S-L-K			
Packing style:	Bulk	Rack	Reloading stack			
Volume Range:	200XL			Variable	Fixed	Multi
	Catalog #					
Standard	AT-200XL-B	AT-200XL-R	AT-200XL-K	2-20 μ l 5-50 μ l 10-100 μ l 20-200 μ l 30-300 μ l	5 μ l, 10 μ l, 20 μ l, 25 μ l, 50 μ l, 100 μ l, 200 μ l	5- 50 μ l 10-100 μ l 20-200 μ l 30-300 μ l
Sterile	AT-200XL-S-B	AT-200XL-S-R	AT-200XL-S-K			
Sterile, Low retention	AT-200XL-S-L-B	AT-200XL-S-L-R	AT-200XL-S-L-K			
Packing style:	Bulk	Rack	Reloading stack			
Volume Range:	1000 μ l			Variable	Fixed	
	Catalog #					
Standard	AT-1000-B	AT-1000-R	AT-1000-K	100-1000 μ l	250 μ l 500 μ l 1000 μ l	
Sterile	AT-1000-S-B	AT-1000-S-R	AT-1000-S-K			
Sterile, Low retention	AT-1000-S-L-B	AT-1000-S-L-R	AT-1000-S-L-K			
Packing style:	Bulk	Rack	Reloading stack			
Volume Range:	1000XL			Variable	Fixed	
	Catalog #					
Standard	AT-1000XL-B	AT-1000XL-R	AT-1000XL-K	100-1000 μ l	250 μ l 500 μ l 1000 μ l	
Sterile	AT-1000XL-S-B	AT-1000XL-S-R	AT-1000XL-S-K			
Sterile, Low retention	AT-1000XL-S-L-B	AT-1000XL-S-L-R	AT-1000XL-S-L-K			
Packing style:	Bulk					
Volume Range:	5000 μ l			Variable		
	Catalog #					
Standard	AT-5000-B			500-5000 μ l		
Sterile	AT-5000B-S-W					
Packing style:	Bulk					
Volume Range:	10000 μ l			Variable	Fixed	Multi
	Catalog #					
Standard	AT-10000-B			1000-10000 μ l		
Sterile	AT-10000B-S-W					

Filter tips				Compatible Volumes		
Packing style:	Bulk	Rack	Reloading stack			
Volume Range:	10 μ l			Variable	Fixed	Multi
	Catalog #					

Standard	AT-10-F-B	AT-10-F-R	AT-10-F-K	0.1 - 2.5µl 0.5-10 µl
Sterile	AT-10-S-F-B	AT-10-S-F-R	AT-10-S-F-K	
Sterile, Low retention	AT-10-S-L-F-B	AT-10-S-L-F-R	AT-10-S-L-F-K	

2.5 µl

0.5 - 10 µl

Packing style:	Bulk	Rack	Reloading stack	
Volume Range:	10XL			Variable
	Catalog #			

Fixed

Multi

Standard	AT-10XL-F-B	AT-10XL-F-R	AT-10XL-F-K	0.1 - 2.5µl 0.5-10 µl
Sterile	AT-10XL-S-F-B	AT-10XL-S-F-R	AT-10XL-S-F-K	
Sterile, Low retention	AT-10XL-S-L-F-B	AT-10XL-S-L-F-R	AT-10XL-S-L-F-K	

2.5 µl

0.5 - 10 µl 2-20 µl

Packing style:	Bulk	Rack	Reloading stack	
Volume Range:	10XN			Variable
	Catalog #			

Fixed

Multi

Standard	AT-10XN-F-B	AT-10XN-F-R	AT-10XN-F-K	0.1 - 2.5µl 0.5-10 µl
Sterile	AT-10XN-S-F-B	AT-10XN-S-F-R	AT-10XN-S-F-K	
Sterile, Low retention	AT-10XN-S-L-F-B	AT-10XN-S-L-F-R	AT-10XN-S-L-F-K	

2.5 µl

0.5 - 10 µl

Packing style:	Bulk	Rack	Reloading stack	
Volume Range:	20µl			Variable
	Catalog #			

Fixed

Multi

Standard	AT-20-F-B	AT-20-F-R	AT-20-F-K	2-20 µl
Sterile	AT-20-S-F-B	AT-20-S-F-R	AT-20-S-F-K	
Sterile, Low retention	AT-20-S-L-F-B	AT-20-S-L-F-R	AT-20-S-L-F-K	

5 µl 10 µl 20 µl

--

Packing style:	Bulk	Rack	Reloading stack	
Volume Range:	50µl			Variable
	Catalog #			

Fixed

Multi

Standard	AT-50-F-B	AT-50-F-R	AT-50-F-K	5-50 µl
Sterile	AT-50-S-F-B	AT-50-S-F-R	AT-50-S-F-K	
Sterile, Low retention	AT-50-S-L-F-B	AT-50-S-L-F-R	AT-50-S-L-F-K	

25 µl 50 µl

5-50 µl

Packing style:	Bulk	Rack	Reloading stack	
Volume Range:	100µl			Variable
	Catalog #			

Fixed

Multi

Standard	AT-100-F-B	AT-100-F-R	AT-100-F-K	10-100 µl
Sterile	AT-100-S-F-B	AT-100-S-F-R	AT-100-S-F-K	
Sterile, Low retention	AT-100-S-L-F-B	AT-100-S-L-F-R	AT-100-S-L-F-K	

100 µl

10 - 100 µl

Packing style:	Bulk	Rack	Reloading stack	
Volume Range:	200µl			Variable
	Catalog #			

Fixed

Multi

Standard	AT-200-F-B	AT-200-F-R	AT-200-F-K	20-200 µl
Sterile	AT-200-S-F-B	AT-200-S-F-R	AT-200-S-F-K	
Sterile, Low retention	AT-200-S-L-F-B	AT-200-S-L-F-R	AT-200-S-L-F-K	

200 µl

20 - 200 µl

Packing style:	Bulk	Rack	Reloading stack	
Volume Range:	200XL			Variable
	Catalog #			

Fixed

Multi

Standard	AT-200XL-F-B	AT-200XL-F-R	AT-200XL-F-K	20-200 µl 30-300µl
Sterile	AT-200XL-S-F-B	AT-200XL-S-F-R	AT-200XL-S-F-K	
Sterile, Low retention	AT-200XL-S-L-F-B	AT-200XL-S-L-F-R	AT-200XL-S-L-F-K	

200 µl

20 - 200 µl 30-300 µl

Packing style:	Bulk	Rack	Reloading stack	
Volume Range:	1000µl			Variable
	Catalog #			

Fixed

Standard	AT-1000-F-B	AT-1000-F-R	AT-1000-F-K	100-1000 µl
Sterile	AT-1000-S-F-B	AT-1000-S-F-R	AT-1000-S-F-K	
Low retention	AT-1000-L-F-B	AT-1000-L-F-R	AT-1000-L-F-K	
Sterile, Low retention	AT-1000-S-L-F-B	AT-1000-S-L-F-R	AT-1000-S-L-F-K	

250 µl 500 µl 1000 µl

Packing style:	Bulk	Rack	Reloading stack	
Volume Range:	1000XL			Variable
	Catalog #			

Fixed

Standard	AT-1000XL-F-B	AT-1000XL-F-R	AT-1000XL-F-K	100-1000 µl
Sterile	AT-1000XL-S-F-B	AT-1000XL-S-F-R	AT-1000XL-S-F-K	
Sterile, Low retention	AT-1000XL-S-L-F-B	AT-1000XL-S-L-F-R	AT-1000XL-S-L-F-K	

250 µl 500 µl 1000 µl

Packing style:	Bulk			
Volume Range:	5000µl			Variable
	Catalog #			

Packing style:	Bulk			
Volume Range:	10000µl			Variable

	Catalog #			
Sterile	AT-1000B-S-F-W			